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
Irritation Assay System (IAS) Ocular and Dermal Irritation

Test Execution and Experimental Protocols

Contents of the Kits

Scientific Background

- The Irritation Assay Systems are in vitro 100% animal free tests that: **Detect**, **Predict** and **Rank** the ocular and dermal irritation potential of several type of materials and substances.


- Cosmetic products
- Household products
- Pharmaceuticals
- Chemical Substances
- Raw Materials

Scientific Background


- The IAS mimic the biochemical phenomena of the ocular and dermal irritation as they would occur in vivo

- The corneal irritancy of a chemical is known to be related to the propensity of said chemical to promote denaturation and disruption of corneal proteins


Scientific Background

Dermal Assay


- Chemicals that cause dermal irritation are known to induce alterations in the structure of keratin, collagen and other dermal proteins.

Scientific Background

- In order to mimic these biochemical phenomena the methods developed make use of 2 important components


Membrane Discs

Proprietary Protein Reagent


Scientific Background – Membrane Discs

- The main difference in the membrane discs between dermal and ocular kits is that the dermal membrane discs have been coated with keratin and collagen in order to more closely mimic the barrier properties of the dermis.
- Additionally a red dye is bound to these proteins during manufacturing, to indicate disruption of the keratin/collagen biobarrier when testing is performed

Scientific Background – Proprietary Protein Reagent

Dermal Assay


Ocular Assay

- The Dermal proprietary protein reagent contains additional globulins with respect to the Ocular one.


- The major constituent is an oligomeric protein consisting of 12 subunits. With a molecular weight of about 320kD, lipids and low molecular weight components

Scientific Background

- When the proprietary protein reagent is rehydrated with a buffering salt solution, proteins, glycoproteins, lipids and other low molecular weight components combine with additional reagent constituents and form an ordered macromolecular matrix.
- Mimicking this way the structure of the Dermis or the Cornea


Reconstituted Protein Solution


Ordered Macromolecular Matrix

Scientific Background

- The application of chemical irritants to the protein reagent solution will lead to protein denaturation. This change in conformation will disrupt the highly ordered macromolecular matrix which will gradually form minute insoluble particles making the clear protein solution cloudy

Ordered Macromolecular Matrix
Clear Solution


Disordered Macromolecular
Matrix
Cloudy Solution

Scientific Background


Ocular Assay

Scientific Background


Dermal Assay


Scientific Background / Quantification

- These phenomena can be quantified by measuring the change in light scattering that occurs as the protein matrix becomes disrupted and the turbidity of the solution increases.
- The light scattering or optical density can be measured by using a plate reader, fitted with the Irritection Assay System software


450nm

405nm


Scientific Background


- Calibrator substances of known irritancy potential are measured together with the test samples.
- This allows the construction of a standard curve that directly relates the Optical Density measurements to an Irritation test score determined by in vivo studies

- The comparison of OD produced by the test sample to the standard curve permits the calculation of an irritancy score that has been shown to be directly related to the potential irritancy of the test material

Scientific Background / Typical Chemical Irritants


Scientific Background / Surfactant Materials

Scientific Background / Irritancy Score

Human Irritancy Equivalent (HIE) Score	Predicted Dermal Irritancy Classification
0.00 - 0.90	Non-Irritant
0.90 - 1.20	Non-Irritant/ Irritant
1.20 - 5.00	Irritant

Dermal Assay


Ocular Assay

Irritation Draize Equivalent (IDE)	Predicted Ocular Irritancy Classification
0.0 – 12.5	Minimal Irritant
>12.5 – 30	Mild Irritant
>30 – 51	Moderate Irritant
> 51	Severe Irritant

Test Execution and Experimental Protocols

Contents of the Kits

Test Execution and Experimental Protocols

Contents of the Kit

Proprietary Reagent Powder	When hydrated forms a solution containing an ordered macromolecular matrix
Hydrating solution	Employed to rehydrate the Reagent Powder and facilitate formation of the ordered protein matrix
Blanking Buffer	Employed as control solution which accounts for the test sample background contribution to the assay
Activator_A	Lowers the pH of the reagent solution the appropriate level to initiate formation of the ordered macromolecular matrix
Calibrators: Cal0, Cal1, Cal2, Cal 3	Known irritants that are employed in each assay to provide standardization and determination of irritancy scores
Quality Controls: QC1 and QC2	Known irritants that are employed in each assay as quality assurance controls to ensure proper performance of the assay
Inhibition Check	Strong irritant substance that is employed as a positive control to check for false negative results at completion of assay
Membrane Discs	Semi permeable membranes that facilitate controlled delivery of the test sample to protein reagent


Kit Formats

Both Ocular and Dermal Irritation kits come in a variety of sizes that can accommodate testing of multiple substances


- 1 sample kit
- 2 samples kit
- 3 samples kit
- 4 samples kit

Test Execution / 1. Preparation of Test Substance


Test Execution / 1. Preparation of Test Substance

Test Execution / 1. Preparation of Test Substance

Bubbles Test

- Solution at 10% of your test substance
- Vortex for 10 seconds
- Allow it to stand for 5 minute
- Examine the sample to see if there is persistent layer of bubbles

Test Execution / 1. Preparation of Test Substance

Non Surfactant Chemical —————> Volume Dependent dose response protocol

Surfactant Chemical —————> Concentration Dependent dose response protocol

Test Execution / 1. Preparation of Test Substance

IAS	Protocol	Doses applied
Dermal	Volume Dependent	25, 50, 75, 100, 125 μL or mg*
Dermal	Concentration Dependent	1, 5, 10, 25, 50%
Ocular	Volume Dependent	25, 50, 75, 100, 125 μL or mg*
Ocular	Concentration Dependent	0,3125, 0,625, 1,25, 2,5, 5%

*In cases of solids or waxy liquids weigh in mg rather than μL

Ocular Irritation / OECD 496 / DB/ALM-157

- On the 24th of October 2019, Ocular Irritation received OECD approval with the protocol number OECD 496

OECD/OCDE

496

Adopted:
24 October 2019

OECD GUIDELINE FOR TESTING OF CHEMICALS

In vitro Macromolecular Test Method for Identifying i) Chemicals Inducing Serious Eye Damage and ii) Chemicals Not Requiring Classification for Eye Irritation or Serious Eye Damage

Test Execution / 2. Reagent Preparation

- Place Hydrating Solution in 25°C incubator and remove membrane discs plate from refrigerator. Allow them to come to room temperature. Do this procedure 2 hours before commencing the test.
- **Rehydration:** Pour all of the hydrating solution into the reagent powder and gently swirl. Let the dissolved reagent stand at RT for about 10 minutes
- **Filtration:** Pour all of the dissolved reagent into a funnel using the filter paper provided within the kits.

Test Execution / 2. Reagent Preparation

- Record the pH of the dissolved and filtered reagent verifying that it falls within the specifications provided in the Range Specification Data Sheet

RANGE SPECIFICATION DATA SHEET

Lot Number: IO 090214

Expiration Date: September 2016


Notice Date: October 2014

pH Before Activation:	7.81 - 8.19
pH After Activation:	5.70 - 5.88
Activator/ Filtered Reagent:	2400 µL/ 40 mL
Activator/ Blanking Buffer:	1800 µL/ 30 mL

- Activation: add the activator reagent to the dissolved and filtered reagent and record the pH once more. Make sure it falls within the specifications of the Range Specification Data Sheet. Add the activator to the Blanking Buffer as well

Test Execution / 3. Test Material Exposure Procedure

- Follow the schemes provided in the instructions manual about how to layout samples, calibrators and quality controls on the 24 well plate.

Calibrators


Quality Controls

Note that for DB-ALM protocol surfactant, the samples are directly pipetted directly in the wells below the membrane discs


Test Execution / 4. Incubation and After

- Cover the 24 well plate with parafilm and incubate for 24 hours.
- After incubation transfer the 24 well assay plate to a 96 well reading plate according to the diagram of the manual
- Place into the plate reader and set up the software for reading.

Irritection Software

- Serves as a user interface to the plate reader
- Processing of spectrophotometric readings
- Compares the increase in OD produced by the samples to the standard curve
- Calculates the “irritancy score” for each tested sample
- Provides graphical representation of the results

ASSAY REPORT - ORIGINAL

Sample Description :		Date :	06/04/15
Sample Number :	15_IT021	Time :	10:53:56
Product Type :	COSMETICO	Technician Name :	LOSINI
Assay Method :	Dermal	Kit Lot Number :	ID120213
Protocol :	Cosmetic_ID120213_04_6	Reagent temperature :	0.0
Incubation Time :	24.0 hours	Reagent pH Before Activation :	10.01
Plate Layout :	5 Samples*4 Volumes	Reagent pH After Activation :	8.30
Instrument Type :	Tecan (auto select)	Sample pH :	
Wavelength :	450nm	Assay Number :	3
Comment :		Assay Qualification :	Qualified

Sample Results:

Dose	Sample OD	Blank OD	Net OD	Irritancy Score	Irritancy Classification	Qualification
50 ul	471	9	462	2.03	Irritant	Qualified
75 ul	565	6	559	2.40	Irritant	Qualified
100 ul	628	6	622	2.64	Irritant	Qualified
125 ul	719	9	710	2.98	Irritant	Qualified

Calibrator Values:

Designation	OD	Irritancy Score	Range Limit (OD)	Qualification
Cal 0	170	0.00	0 - 200	Range qualified
Cal 1	196	1.00	104 - 260	Range qualified
Cal 2	454	2.00	330 - 630	Range qualified
Cal 3	977	4.00	810 - 1430	Range qualified

Quality Control Values:

Designation	OD	Irritancy Score	Range Limit (Score)	Qualification
QC 1	147	0.75	0.11 - 0.95	Range qualified
QC 2	723	3.03	0.94 - 3.80	Range qualified

* Mean value from assay data history

** Mean value from protocol defaults or adjusted value due to calibrator zero substitution

[] Value before substitution

Irritection
Software /
Qualified
Assay

Results Interpretation and Data Analysis

Indicates that values obey expected behaviour of typical dose response curve

Dose	Sample OD	Blank OD	Net OD	Irritancy Score	Irritancy Classification	Qualification
25 mg	124	-4	128	7.4	Minimal	Qualified
50 mg	114	-1	115	6.7	Minimal	Qualified
75 mg	140	-2	142	8.2	Minimal	Qualified
100 mg	151	-4	155	9.0	Minimal	Qualified
125 mg	134	-5	139	8.0	Minimal	Qualified

The highest irritancy score calculated by the software is the irritancy of the test sample and it is called Maximum Qualified Score (MQS)
 It is the value that correlates the most with the in vivo irritancy property of the test material


Results Interpretation and Data Analysis / Reporting

Sample	Protocol	Method	Max IDE Score	Predicted Dermal Irritancy Classification
Campione 1	Cosmetic / Volume dependent dose response curve	Ocular Irritection	9.0	Minimal Irritant

Interpretation of Ocular Irritation[®] Results According to Protocol

Protocol DIE	DB-ALM 157	OECD 496	GHS Classification
0 – 12,5	Minimal	No Category	Non-Irritant
> 12,5 – 30	Mild	No Prediction	Irritant
> 30 – 51	Moderate	Category 1	
> 51	Severe		

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